

# Meeting Minutes - Cut Tag Project Report

---

- **Date:** 2025-08-30

## 1. Nucleosome Architecture and Histone Organization

- **Key Points**

- The nucleosome core particle consists of ~147 bp of DNA wrapped ~1.7 turns around a histone octamer.
- **Histone spatial arrangement:**
  - Two copies each of H2A, H2B, H3, and H4 form the octamer.
  - The central scaffold is the (H3–H4)<sub>2</sub> tetramer; two H2A–H2B dimers dock on either side to complete the octamer.
  - DNA engages the histone octamer via the acidic patch (notably on H2A–H2B) and multiple basic histone tails.
- **Histone H1 (linker histone):**
  - Binds near the nucleosome dyad and the entering/exiting linker DNA.
  - **Stabilizing roles:**
    - Enhances DNA wrapping and reduces nucleosome breathing.

## 2. CUT&Tag vs ChIP-seq: Methods and Practical Considerations

- **Method Overviews**

- **ChIP-seq** (Chromatin Immunoprecipitation followed by sequencing):
  - Workflow: Crosslink (often formaldehyde) → chromatin fragmentation (sonication or MNase) → immunoprecipitation with antibody → reverse crosslinks → library prep → sequencing.
- **CUT&Tag** (Cleavage Under Targets and Tagmentation):
  - Workflow: Lightly permeabilize nuclei/cells → primary antibody binds target → secondary antibody recruits Protein A/G–Tn5 transposase loaded with adapters → tethered Tn5 cleaves and tags DNA at antibody-bound sites → direct PCR amplification → sequencing.

- **Comparative Discussion**

- **Resolution and signal-to-noise:**
  - **CUT&Tag:** High resolution and excellent signal-to-noise due to tethered cleavage and minimal background.
  - **ChIP-seq:** Broader fragments; background from nonspecific fragmentation and immunoprecipitation.
- **DNA fragment characteristics:**
  - **CUT&Tag:** Generates DNA fragments of highly consistent and appropriate lengths, which significantly improves the efficiency and accuracy of *de novo* DNA sequencing.
  - **ChIP-seq:** Produces a broader range of fragment sizes due to nonspecific sonication, complicating subsequent sequencing and assembly.
- **Workflow complexity and turnaround:**
  - **CUT&Tag:** Fewer steps; direct library generation via Tn5 reduces hands-on time.

- **ChIP-seq:** More steps; longer turnaround.
- **Compatibility:**
  - **CUT&Tag:** Works well for histone marks and many transcription factors; dependent on antibody quality and accessibility.

### 3. Nucleic Acid Chemistry, Directionality, and Nomenclature

- **DNA structure and numbering**

- Backbone: Repeating deoxyribose and phosphate groups; bases attached at C1'.
- Sugar numbering: C1' (base), C2' (H in DNA; OH in RNA), C3' (OH), C4', C5' (phosphoryl linkage).
- Strands are antiparallel with polarity (5' → 3').

- **Phosphates on nucleotide building blocks**

- Clarification from our discussion:
  - The phosphodiester backbone between nucleotides contains one phosphate per linkage.
  - The 5' end of a nucleic acid chain typically bears a monophosphate (5'-P) or hydroxyl, depending on preparation.
  - The incoming building blocks for polymerization are nucleoside triphosphates:
    - For DNA synthesis: dNTPs (dATP, dCTP, dGTP, dTTP).
    - For RNA synthesis: NTPs (ATP, CTP, GTP, UTP).
  - These carry a triphosphate on the 5' carbon; during incorporation, pyrophosphate (PPi) is released, leaving a single phosphate in the backbone.
  - Conclusion recorded: Each incorporated nucleotide contributes one phosphate in the backbone; the free nucleotide substrates are triphosphates prior to incorporation.

- **Abbreviations**

- dNTP: Deoxyribonucleoside triphosphate (dATP, dCTP, dGTP, dTTP).
- ATP/GTP/CTP/TTP: Ribonucleoside triphosphates for RNA (note: UTP replaces TTP in RNA; TTP is typically used in DNA contexts as dTTP).

- **Template vs non-template strands**

- Template strand: The DNA strand read by polymerases to synthesize a complementary strand.
- Non-template (coding) strand: Same sequence as the RNA transcript (except T ↔ U), not directly read during transcription.

- **Directional conventions in core processes**

- DNA replication: DNA polymerases extend DNA only (5' → 3'); leading strand continuous, lagging strand discontinuous (Okazaki fragments).
- Transcription: RNA polymerase synthesizes RNA (5' → 3') while reading the DNA template (3' → 5').
- Translation: Ribosomes read mRNA (5' → 3'); polypeptide synthesized from N-terminus to C-terminus.

### 4. DNA Sequencing Quality Metrics and GC Content

- **Q-scores (Phred quality)**

- Definition:  $(Q = -10 \log_{10}(P_{\text{error}}))$ .
- Common thresholds:
  - Q20:  $(P_{\text{error}} = 1\%)$  (accuracy ~99%).
  - Q30:  $(P_{\text{error}} = 0.1\%)$  (accuracy ~99.9%).
- Usage:
  - Reported as per-base scores and aggregate metrics (e.g.,  $\% \geq Q30$ ).
  - Higher  $\% \geq Q30$  indicates better overall run quality; filters for downstream analyses.

- **GC content**

- Meaning: Fraction of G and C bases in reads or genome regions.
- Why we check it:
  - Detects library prep or sequencing biases (e.g., drop-outs at extreme GC).
  - Flags contamination or sample swaps if observed GC deviates from expected.
  - Informs normalization and peak-calling (AT/GC bias corrections) and guides PCR conditions.